

## Calciums on Beckman Instruments

It has been suggested to us that one of the reasons why calcium is apparently not performed as well as other chemistries, is due to issues surrounding registration. In order to investigate this, we would like all Beckman instrument users, to complete the form below. The information will allow us to investigate the statistics around the various ways in which calcium is handled by Beckman instruments.

We will forward the outcome of this mini-survey to all Beckman users, as soon as we have completed the stats investigation.

Lab QA No: \_\_\_\_\_ Lab Name: \_\_\_\_\_

Lab email address: \_\_\_\_\_

Lab contact person: \_\_\_\_\_

### *Which Beckman group do you belong to?*

- Group 1: Indirect ISE  
- CX 3 Delta, CX 5 Delta, CX 9 (ALX 9), CX 5/9 Pro
- Group 2: Arsenazo (Cup Chemistry)  
- CX 3, CX 7
- Group 3: Arsenazo (Cartridge Chemistry)  
- CX 4/5/7/9 CE, CX Delta, CX 9 (ALX 9), CX 5/9 Pro
- Group 4: Indirect ISE  
LX 20, DxC 600, DxC 800

Dr Jim McCulloch  
July 2006

S A N A S



PROFICIENCY TESTING Accredited to ISO Guide 43 and ILAC G13